

| | | |
|--------------------------------|---|--|
| Product Name: | Examination Gloves | |
| Product Application: | NEST gloves can be used in household, industrial, medical, beauty and other industries, and it is an essential hand protection product. | |
| Product Materials: | Latex / Nitrile + colloidal oatmeal | |
| product Specifications: | Temperature range: stored at room temperature Shelf life: 3 Years after date of production (ensure package is in good) | |
| Sterilization: | No | |
| BSE/TSE: | These products are deemed animal free. | |
| Performance Testing: | Each manufacturing lot is sampled and tested in accordance with standard operating procedures. Appearance inspection: qualified Tensile Strength Test: qualified Leak test: qualified Packaging inspection: qualified | |
| Features: | <ol style="list-style-type: none"> 1. Single use only 2. Made in Malaysia 3. Ambidextrous: Can be worn left and right, regardless of left or right hand 4. Glove mouth rolled edge design, easy to wear, prevent rollback 5. Textured outer surface of the glove increases friction and facilitates gripping 6. Powder-Free 7. Non-Sterile 8. Boxed for easy extraction. 9. Available in three materials <ol style="list-style-type: none"> (1) Nitrile Examination Glove: <ul style="list-style-type: none"> • Tensile Strength: Minimum 18.0 MPa (Before Aging) • Palm thickness: 0.07 mm (2) Latex Examination Glove: <ul style="list-style-type: none"> • Natural rubber latex and compounding materials • Tensile Strength: Minimum 18.0 MPa (Before Aging) • Palm thickness: 0.11 mm (3) Colloidal Oatmeal Coated Nitrile Glove Colloidal Oatmeal Coated Latex Glove <ul style="list-style-type: none"> • Oats Extractions <p>Oats extractions is a patented coating recognized by the FDA as a skin protectant, which have been proven by clinical studies could reduce and control rash and dry skin symptoms. Clinical studies have shown that colloidal oatmeal could reduce and control rash, dry skin.</p> | |

Product Range

Nitrile Examination Glove

| Cat. No. | Size | Palm Width (mm) ± 4 | Length | /Inner Box | /Case |
|----------|------|-------------------------|---------------------|------------|-------|
| 902001 | XS | 76 | Min 230mm (ASTM) | 100 | 1000 |
| 902011 | S | 86 | | 100 | 1000 |
| 902021 | M | 98 | | 100 | 1000 |
| 902031 | L | 107 | | 100 | 1000 |
| 902041 | XL | 115 | | 90 | 900 |

Latex Examination Glove

| Cat. No. | Size | Palm Width (mm) ± 4 | Length | /Inner Box | /Case |
|----------|------|-------------------------|---------------------|------------|-------|
| 901001 | XS | 76 | Min 230mm (ASTM) | 100 | 1000 |
| 901011 | S | 86 | | 100 | 1000 |
| 901021 | M | 98 | | 100 | 1000 |
| 901031 | L | 107 | | 100 | 1000 |
| 901041 | XL | 115 | | 90 | 900 |

Colloidal Oatmeal Coated Nitrile Glove

| Cat. No. | Size | Palm Width (mm) ± 4 | Length | /Inner Box | /Case |
|----------|------|-------------------------|---------------------|------------|-------|
| 903011 | S | 84 | Min 230mm (ASTM) | 100 | 1000 |
| 903021 | M | 94 | | 100 | 1000 |
| 903031 | L | 105 | | 100 | 1000 |

Colloidal Oatmeal Coated Latex Glove

| Cat. No. | Size | Palm Width (mm) ± 4 | Length | /Inner Box | /Case |
|----------|------|-------------------------|---------------------|------------|-------|
| 904011 | S | 84 | Min 230mm (ASTM) | 100 | 1000 |
| 904021 | M | 94 | | 100 | 1000 |
| 904031 | L | 105 | | 100 | 1000 |

Product Picture

Nitrile Examination Glove



Due to chromatic aberration of the display screen, there may be chromatic aberration between the picture and the real object.

Latex Examination Glove



The picture is only for reference, please make the object as the standard.