

Certificate of Analysis

Product Name: Polyclonal Anti-wfikkn1 Antibody

Catalog Number: DZ11633

Description

Lot No.	19DZ47F03
Size	100µl/vial
Form	Liquid
Concentration	1mg/ml
Ig type	Rabbit IgG
Purification	Immunogen affinity purified.
Contents	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .

Data of Indirect ELISA

Coating concentration of antigen	0.1µg/ml	0.1µg/ml	0.1µg/ml	0.1µg/ml
Working concentration of anti- wfikkn1 antibody	1µg/ml	0.1µg/ml	0.01µg/ml	0.001µg/ml
Working concentration of HRP conjugated goat anti-rabbit IgG (Catalog # BA1054)	0.2µg/ml	0.2µg/ml	0.2µg/ml	0.2µg/ml
O.D. values (450nm)	1.856	0.603	0.137	0.050
Antibody Status	Pass			


Data Evaluation Criteria

Coating concentration of antigen	0.1µg/ml	0.1µg/ml
Working concentration of the primary antibody	1µg/ml	1µg/ml
Working concentration of HRP conjugated goat anti-rabbit IgG	0.2µg/ml	0.2µg/ml
O.D. values (450nm)	≥ 1.0	<1.0
Antibody's affinity for the antigen	Strong	Weak
Antibody Status	Pass	Fail

Storage

At -20°C for one year, at 4°C for one month. Avoid repeated freezing and thawing.

Approved by:



Ms. Yangjing Yue

Chief technician

Quality Control System of ELISA

Date: June 3, 2019