


AMPK Signaling Pathway



Key Products for Energy Metabolism

Visit www.bosterbio.com for full listing of antibodies, related kits and reagents.

ANTIBODIES		APPLICATION	REACTIVITY
PB9418	Anti-Adiponectin Receptor 1/ADIPOR1 Picoband™ Antibody	WB	Human, Rat
M00024-1	Anti-AKT1 Monoclonal Antibody	FC, IP, IF, IHC, WB	Human, Mouse, Rat
PB9043	Anti-AKT2 Picoband™ Antibody	IHC, WB	Human, Mouse, Rat
PB9752	Anti-alpha 1d Adrenergic Receptor/ADRA1A Picoband™ Antibody	WB	Human, Mouse, Rat
A00994-3	Anti-AMPK alpha 1/PRKAA1 Picoband™ Antibody	ELISA, WB	Human, Mouse, Rat
PB9738	Anti-AMPK beta 1/PRKAB1 Picoband™ Antibody	WB	Human
M05077	Anti-AMPK beta 2 Picoband™ Antibody (monoclonal, 6G1)	FC, IHC, ICC, WB	Human
PB9739	Anti-AMPK beta 2/PRKAB2 Picoband™ Antibody	IHC, WB	Human, Mouse, Rat
PB9481	Anti-ATG14L Picoband™ Antibody	FC, IHC, ICC, WB	Human, Rat
PB9076	Anti-Becn1 1/BECN1 Picoband™ Antibody	IHC, WB	Human, Mouse, Rat
PA1520	Anti-CaMKK/CAMKK1 Antibody	IHC, WB	Human, Mouse, Rat
PA1497	Anti-CAMKK2 Antibody	IHC, ICC, WB	Human, Mouse, Rat
PB9491	Anti-CPT1B Picoband™ Antibody	IHC, WB	Human, Mouse, Rat
PB9100	Anti-CREB/CREB1 Picoband™ Antibody	WB	Human, Rat
PB9540	Anti-Cryptochrome I/CRY1 Picoband™ Antibody	IHC, WB	Human, Mouse, Rat
PB9865	Anti-Fatty Acid Synthase/FASN Picoband™ Antibody	IHC, ICC, WB	Human
PB9195	Anti-FOXO1A Picoband™ Antibody	WB	Human, Rat
PB9196	Anti-FOXO3A Picoband™ Antibody	IHC, WB	Human, Rat
PB9109	Anti-Glucose Transporter GLUT4/SLC2A4 Picoband™ Antibody	IHC, WB	Human, Mouse, Rat
PA2012	Anti-Hamartin/TSC1 Antibody	IHC, WB	Human, Mouse, Rat
A00971-1	Anti-HDAC4 Picoband™ Antibody	WB	Human, Mouse, Rat
A01230-5	Anti-HDAC5 Picoband™ Antibody	ELISA, IHC-P, WB	Human, Mouse, Rat
PB9215	Anti-HNF-4-alpha/HNF4A Picoband™ Antibody	FC, IF, IHC-P, IHC-F, ICC, WB	Human, Mouse, Rat
PA1620	Anti-Insulin Receptor/INSR Antibody	WB	Human, Mouse, Rat
PB9178	Anti-KAT3B/p300/EP300 Picoband™ Antibody	IHC, WB	Human, Mouse, Rat
PB9480	Anti-KIAA0652/ATG13 Picoband™ Antibody	WB	Human, Mouse
A00350	Anti-Leptin Receptor/LEPR Picoband™ Antibody	ELISA, IHC	Human, Mouse, Rat
PA1355	Anti-LKB1/STK11 Antibody	IHC, WB	Human, Mouse, Rat
PB9737	Anti-Myosin Phosphatase/PPP1R12A Picoband™ Antibody	FC, IHC, ICC, WB	Human, Mouse
P00024-5	Anti-Phospho-Akt (Ser473) Rabbit Monoclonal Antibody	IF, IHC, ICC, WB	Human, Mouse, Rat
PB9346	Anti-Phospholipase C beta 1/PLCB1 Picoband™ Antibody	WB	Human, Mouse, Rat
A01091-1	Anti-PI3 Kinase p110 beta/PIK3CB Picoband™ Antibody	IHC, WB	Human, Mouse, Rat
A06142-1	Anti-RSK3/RPS6KA2 Picoband™ Antibody	ELISA, IF, IHC-P, ICC, WB	Human, Mouse, Rat
PA1482	Anti-Tau/MAPT Antibody	WB	Human, Mouse, Rat
A01118-1	Anti-TORC2/CRTC2 Picoband™ Antibody	WB	Human, Mouse, Rat
PB9121	Anti-Tuberin/TSC2 Picoband™ Antibody	IHC, WB	Human, Mouse, Rat


APPLICATIONS: WB Western Blotting / IP Immunoprecipitation / IHC Immunohistochemistry / IF Immunofluorescence / FC Flow Cytometry / ICC Immunocytochemistry
ChIP Chromatin Immunoprecipitation / -P Paraffin / -F Frozen

 **25**
Years

 **23,000+**
Publications

 **4.8/5**
Rating on biocompare.com

 **16,000+**
Antibodies

 **1000+**
ELISA Kits

BOSTER
antibody and ELISA experts



For more details, visit www.bosterbio.com
3942 B Valley Ave, Pleasanton, CA, 94566 | Toll free: 888-466-3604
Email: support@bosterbio.com