


Hypoxia Signaling Pathway



Key Products for Cellular Metabolism

Visit www.bosterbio.com for full listing of antibodies, related kits and reagents.

ANTIBODIES		APPLICATION	REACTIVITY
A00024	Anti-AKT1 Picoband™ Antibody	WB	Mouse, Rat
PB9043	Anti-AKT2 Picoband™ Antibody	IHC, WB	Human, Mouse, Rat
A00994-3	Anti-AMPK alpha 1/PRKAA1 Picoband™ Antibody	ELISA, WB	Human, Mouse, Rat
PB9738	Anti-AMPK beta 1/PRKAB1 Picoband™ Antibody	WB	Human
PB9739	Anti-AMPK beta 2/PRKAB2 Picoband™ Antibody	IHC, WB	Human, Mouse, Rat
PB9164	Anti-B Raf/BRAF Picoband™ Antibody	WB	Human, Mouse, Rat
M00004-2	Anti-beta Catenin Picoband™ Antibody (monoclonal, 1F6)	FC, IHC, ICC, WB	Human, Mouse, Rat
A00004	Anti-beta Catenin/CTNNB1 Picoband™ Antibody	ELISA, FC, IF, IHC-P, IHC-F, ICC, WB	Human, Mouse, Rat
PA1343	Anti-ERK1/MAPK3 Antibody	IHC, ICC, WB	Human, Mouse, Rat
PA1310-1	Anti-ERK2/MAPK1 Antibody	WB	Human, Mouse, Rat
PB9129	Anti-HIF-1 beta/ARNT Picoband™ Antibody	WB	Human
PB9253	Anti-HIF-1-alpha/HIF1A Picoband™ Antibody	IHC, WB	Human, Mouse, Rat
PA1129-2	Anti-HIF-2-alpha/EPAS1 Antibody	WB	Human, Mouse, Rat
PA2145	Anti-HIF-3/HIF3A Antibody	WB	Human, Rat
PB9110	Anti-IKK alpha/CHUK Picoband™ Antibody	IHC, WB	Human, Mouse, Rat
PB9221	Anti-IKK beta/IKBKB Picoband™ Antibody	IHC, WB	Human
A00874	Anti-IKK gamma/IKBKG Picoband™ Antibody	WB	Human, Mouse, Rat
A00330	Anti-JAK1 Picoband™ Antibody	WB	Human, Mouse, Rat
A00284-1	Anti-NF-kB p65/RELA Picoband™ Antibody	ELISA, IHC, WB	Human, Mouse, Rat
PB9008	Anti-P53/TP53 Picoband™ Antibody	IHC, WB	Human
A02182-1	Anti-PHD3/EGLN3 Picoband™ Antibody	IHC, WB	Human, Mouse, Rat
MA1055	Anti-Phospho-MAP Kinase, Activated Antibody (Monoclonal, MAPK-YT)	IHC, ICC, WB	Human, Mouse, Rat, Yeast
A00446-1	Anti-Raf1 Picoband™ Antibody	WB	Human, Mouse, Rat
PB9159	Anti-SIRT1 Picoband™ Antibody	WB	Human
PB9357	Anti-SIRT6 Picoband™ Antibody	IHC, WB	Human, Mouse
PA1108	Anti-STAT3 Antibody	IHC, ICC, WB	Human, Mouse, Rat
A05315	Anti-TCR alpha/TRAC Picoband™ Antibody	ELISA, FC, WB	Human, Mouse, Rat
A00429-2	Anti-TLR1 Picoband™ Antibody	ELISA, FC, IHC, ICC, WB	Human, Mouse, Rat
A00131	Anti-TLR2 Picoband™ Antibody	WB	Human, Mouse, Rat
A00197	Anti-TLR3 Picoband™ Antibody	ELISA, WB	Human, Mouse, Rat
PA1484	Anti-TLR4 Antibody	IHC, WB	Human
PB9451	Anti-TLR5 Picoband™ Antibody	FC, IHC, ICC, WB	Human
PB9452	Anti-TLR7 Picoband™ Antibody	FC, IHC, ICC, WB	Human, Mouse, Rat
RP1073	Anti-TLR8 Antibody	WB	Human, Rat
PA1958	Anti-TLR10 Antibody	WB	Human
PB9827	Anti-Transferrin/TF Picoband™ Antibody	IHC, WB	Human, Mouse, Rat


APPLICATIONS: WB Western Blotting / IP Immunoprecipitation / IHC Immunohistochemistry / IF Immunofluorescence / FC Flow Cytometry / ICC Immunocytochemistry
ChIP Chromatin Immunoprecipitation / -P Paraffin / -F Frozen

 **25**
Years

 **23,000+**
Publications

 **4.8/5**
Rating on biocompare.com

 **16,000+**
Antibodies

 **1000+**
ELISA Kits

BOSTER
antibody and ELISA experts



For more details, visit www.bosterbio.com
3942 B Valley Ave, Pleasanton, CA, 94566 | Toll free: 888-466-3604
Email: support@bosterbio.com