


# Insulin Receptor Signaling Pathway



# Key Products for Cellular Metabolism

Visit [www.bosterbio.com](http://www.bosterbio.com) for full listing of antibodies, related kits and reagents.

ANTIBODIES		APPLICATION	REACTIVITY
A00024	Anti-AKT1 Picoband™ Antibody	WB	Mouse, Rat
PB9043	Anti-AKT2 Picoband™ Antibody	IHC, WB	Human, Mouse, Rat
M00725-1	Anti-AKT2 Picoband™ Antibody (monoclonal, 10C6)	ICC, WB	Human
PB9165	Anti-Caveolin-1/CAV1 Picoband™ Antibody	IF, IHC, ICC, WB	Human, Mouse, Rat
PB9166	Anti-Caveolin-2/CAV2 Picoband™ Antibody	IF, IHC, WB	Human, Mouse
A00990-2	Anti-Caveolin-3/CAV3 Picoband™ Antibody	ELISA, IHC-P, WB	Human, Mouse, Rat
PA1800	Anti-CBL Antibody	WB	Human, Mouse, Rat
PA2033	Anti-Flotillin 1/FLOT1 Antibody	WB	Human, Mouse, Rat
PB9623	Anti-Flotillin 2/FLOT2 Picoband™ Antibody	IHC, WB	Human, Mouse, Rat
PB9195	Anti-FOXO1A Picoband™ Antibody	WB	Human, Rat
PB9196	Anti-FOXO3A Picoband™ Antibody	IHC, WB	Human, Rat
A01989-1	Anti-GAB1 Picoband™ Antibody	WB	Human, Mouse, Rat
PB9109	Anti-Glucose Transporter GLUT4/SLC2A4 Picoband™ Antibody	IHC, WB	Human, Mouse, Rat
PA1580-1	Anti-GRB10 Antibody	IHC, WB	Human, Mouse, Rat
PA1579	Anti-GRB2 Antibody	WB	Human, Mouse, Rat
A03152-2	Anti-GSK3 alpha/GSK3A Picoband™ Antibody	WB	Human, Mouse, Rat
PA2012	Anti-Hamartin/TSC1 Antibody	IHC, WB	Human, Mouse, Rat
A00368	Anti-iNOS/NOS2 Picoband™ Antibody	WB	Human
PA1620	Anti-Insulin Receptor/INSR Antibody	WB	Human, Mouse, Rat
PA1892	Anti-JNK1/MAPK8 Antibody	WB	Human, Mouse, Rat
RP1077	Anti-JNK2/MAPK9 Antibody	WB	Human
PA1355	Anti-LKB1/STK11 Antibody	IHC, WB	Human, Mouse, Rat
PA1376	Anti-MEK1/MAP2K1 Antibody	IHC, ICC, WB	Human, Mouse, Rat
A01268-1	Anti-Mitochondrial Pyruvate dehydrogenase kinase 1/PDK1 Picoband™ Antibody	ELISA, IHC, WB	Human, Mouse, Rat
PB9314	Anti-PI 3 Kinase catalytic subunit alpha/PIK3CA Picoband™ Antibody	WB	Human, Rat
PA1410	Anti-PI 3 Kinase p85 alpha/PIK3R1 Antibody	WB	Human, Mouse, Rat
PA1680	Anti-PI 3 Kinase p85 beta/PIK3R2 Antibody	IHC, WB	Human, Mouse, Rat
A01091-1	Anti-PI3 Kinase p110 beta/PIK3CB Picoband™ Antibody	IHC, WB	Human, Mouse, Rat
A01890	Anti-PKC gamma/PRKCG Picoband™ Antibody	IHC, WB	Human, Mouse, Rat
PB9431	Anti-Rac1 Picoband™ Antibody	WB	Human, Mouse, Rat
PA2184	Anti-SGK1 Antibody	IHC, ICC, WB	Human, Mouse, Rat
M00150-2	Anti-SHP2/PTPN11 Picoband™ Antibody (monoclonal, 2E6)	FC, IF, IHC-P, ICC, WB	Human, Mouse, Rat
PA1707	Anti-SOCS3 Antibody	IHC, WB	Human, Mouse
A00294-3	Anti-TNF Receptor I/TNFRSF1A Picoband™ Antibody	ELISA, FC, IHC, ICC, WB	Human, Mouse, Rat
A00294	Anti-TNFRSF1A Picoband™ Antibody	WB	Human, Mouse, Rat
PB9121	Anti-Tuberin/TSC2 Picoband™ Antibody	IHC, WB	Human, Mouse, Rat


APPLICATIONS: WB Western Blotting / IP Immunoprecipitation / IHC Immunohistochemistry / IF Immunofluorescence / FC Flow Cytometry / ICC Immunocytochemistry  
ChIP Chromatin Immunoprecipitation / -P Paraffin / -F Frozen

 **25**  
Years

 **23,000+**  
Publications

 **4.8/5**  
Rating on biocompare.com

 **16,000+**  
Antibodies

 **1000+**  
ELISA Kits

**BOSTER**  
antibody and ELISA experts



For more details, visit [www.bosterbio.com](http://www.bosterbio.com)  
3942 B Valley Ave, Pleasanton, CA, 94566 | Toll free: 888-466-3604  
Email: [support@bosterbio.com](mailto:support@bosterbio.com)