


# Warburg Effect Signaling Pathway



# Key Products for Cellular Metabolism

Visit [www.bosterbio.com](http://www.bosterbio.com) for full listing of antibodies, related kits and reagents.

ANTIBODIES	APPLICATION	REACTIVITY
A00994-3 Anti-AMPK alpha 1/PRKAA1 Picoband™ Antibody	ELISA, WB	Human, Mouse, Rat
PB9738 Anti-AMPK beta 1/PRKAB1 Picoband™ Antibody	WB	Human
M05077 Anti-AMPK beta 2 Picoband™ Antibody (monoclonal, 6G1)	FC, IHC, ICC, WB	Human
PB9739 Anti-AMPK beta 2/PRKAB2 Picoband™ Antibody	IHC, WB	Human, Mouse, Rat
PA1013-1 Anti-Bax Antibody	FC, IHC, ICC, WB	Human, Mouse, Rat
A00183 Anti-Bax Picoband™ Antibody	FC, IHC, ICC, WB	Human, Mouse, Rat
MA1028 Anti-c-Myc Antibody (Monoclonal, 9E10)	IHC, ICC, WB	Human
PB9092 Anti-c-Myc Picoband™ Antibody	IHC-P, WB	Human, Rat
PA1343 Anti-ERK1/MAPK3 Antibody	IHC, ICC, WB	Human, Mouse, Rat
PA1310-1 Anti-ERK2/MAPK1 Antibody	WB	Human, Mouse, Rat
PB9435 Anti-Glucose Transporter GLUT1/SLC2A1 Picoband™ Antibody	WB	Human, Mouse, Rat
PB9436 Anti-Glucose Transporter GLUT2/SLC2A2 Picoband™ Antibody	IHC, WB	Human
PB9109 Anti-Glucose Transporter GLUT4/SLC2A4 Picoband™ Antibody	IHC, WB	Human, Mouse, Rat
PB9960 Anti-Glucose Transporter GLUT5/SLC2A5 Picoband™ Antibody	FC, IF, IHC-P, IHC-F, ICC, WB	Human, Rat
PB9891 Anti-GLUT9/SLC2A9 Picoband™ Antibody	IHC, WB	Human, Mouse, Rat
PA2167 Anti-GLUT12/SLC2A12 Antibody	WB	Human, Mouse, Rat
A01389 Anti-Hexokinase II/HK2 Picoband™ Antibody	WB	Human, Mouse, Rat
PB9253 Anti-HIF-1-alpha/HIF1A Picoband™ Antibody	IHC, WB	Human, Mouse, Rat
PB9602 Anti-IDH2 Picoband™ Antibody	IHC, ICC, WB	Human, Mouse, Rat
PB9601 Anti-Isocitrate dehydrogenase/IDH1 Picoband™ Antibody	FC, IF, IHC-P, IHC-F, ICC, WB	Human, Mouse, Rat
PB10075 Anti-LDHA Picoband™ Antibody	IHC, WB	Human, Mouse, Rat
PA1355 Anti-LKB1/STK11 Antibody	IHC, WB	Human, Mouse, Rat
A01268-1 Anti-Mitochondrial Pyruvate dehydrogenase kinase 1 Picoband™ Antibody	ELISA, IHC, WB	Human, Mouse, Rat
MA1078 Anti-P53 Monoclonal Antibody	IHC, ICC, WB	Human
PB9008 Anti-P53/TP53 Picoband™ Antibody	IHC, WB	Human
PA1994 Anti-PDK2 Antibody	IHC, ICC, WB	Human, Mouse, Rat
A03261 Anti-PDK2 Picoband™ Antibody	ELISA, WB	Human, Mouse, Rat
PB9773 Anti-PDK4 Picoband™ Antibody	WB	Human, Mouse, Rat
PB9314 Anti-PI 3 Kinase catalytic subunit alpha/PIK3CA Picoband™ Antibody	WB	Human, Rat
PA1410 Anti-PI 3 Kinase p85 alpha/PIK3R1 Antibody	WB	Human, Mouse, Rat
PA1680 Anti-PI 3 Kinase p85 beta/PIK3R2 Antibody	IHC, WB	Human, Mouse, Rat
A01091-1 Anti-PI3 Kinase p110 beta/PIK3CB Picoband™ Antibody	IHC, WB	Human, Mouse, Rat
PB9379 Anti-PKM2 Picoband™ Antibody	IHC, WB	Human
M00163-1 Anti-SLC2A1 Picoband™ Antibody (monoclonal, 10C10)	IHC-P, WB	Human
A02256-1 Anti-UCP2 Picoband™ Antibody	WB	Human, Mouse, Rat
PA1780 Anti-VDAC/Porin/VDAC1 Antibody	IHC, ICC, WB	Human, Mouse, Rat


APPLICATIONS: WB Western Blotting / IP Immunoprecipitation / IHC Immunohistochemistry / IF Immunofluorescence / FC Flow Cytometry / ICC Immunocytochemistry  
ChIP Chromatin Immunoprecipitation / -P Paraffin / -F Frozen

 **25**  
Years

 **23,000+**  
Publications

 **4.8/5**  
Rating on *biocompare.com*

 **16,000+**  
Antibodies

 **1000+**  
ELISA Kits

**BOSTER**  
antibody and ELISA experts



For more details, visit [www.bosterbio.com](http://www.bosterbio.com)  
3942 B Valley Ave, Pleasanton, CA, 94566 | Toll free: 888-466-3604  
Email: [support@bosterbio.com](mailto:support@bosterbio.com)