BOSTER BIOLOGICAL TECHNOLOGY
3942 B Valley Ave, Pleasanton, CA 94566
888-466-3604 | support@bosterbio.com | www.bosterbio.com

Anti-TFIIIC102 (E131) GTF3C3 Antibody

Catalog Number: A12003-1

Submit a product review to Biocompare.com





Submit a review of this product to Biocompare.com to receive a $20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Overview


Product Name	Anti-TFIIIC102 (E131) GTF3C3 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-TFIIIC102 (E131) GTF3C3 Antibody catalog # A12003-1. Tested in IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IHC
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9Y5Q9

Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 100-150 of Human TFIIIC102.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	IHC: 1:50-1:200 This protein A purified antibody has been tested for use in ELISA, immunofluorescence microscopy and western blotting. Specific conditions for reactivity should be optimized by the end user. Expect a band approximately 26-31 kDa in size corresponding to human CENP-Q by western blotting in the appropriate cell lysate or extract.

BOSTER BIOLOGICAL TECHNOLOGY
3942 B Valley Ave, Pleasanton, CA 94566
888-466-3604 | support@bosterbio.com | www.bosterbio.com

Anti-TFIIIC102 (E131) GTF3C3 Antibody (A12003-1) Images



Immunohistochemistry (IHC) analyzes of TFIIIC102 (E131) pAb in paraffin-embedded human breast carcinoma tissue at 1:50.

Anti-TFIIIC102 (E131) GTF3C3 Antibody

For Research Use Only. Not for use in diagnostic procedures.