BOSTER BIOLOGICAL TECHNOLOGY
3942 B Valley Ave, Pleasanton, CA 94566
888-466-3604 | support@bosterbio.com | www.bosterbio.com

Tissue Protein Lysate, Colon: right

Catalog number: LS781203

This package insert must be read in its entirety before using this product.

For research use only. Not for use in diagnostic procedures.

Tissue Protein Lysate, Colon: right

BOSTER BIOLOGICAL TECHNOLOGY
3942 B Valley Ave, Pleasanton, CA 94566
888-466-3604 | support@bosterbio.com | www.bosterbio.com

Catalog Number: LS781203, **Storage:** The lysate is shipped with dry ice. Store the tissue lysate at -80°C. These products are stable for one year from date of shipping when stored at -80°C without repeated freeze/thaws.

Overview

Product Name	Tissue Protein Lysate, Colon: right
Sku	LS781203
Size	150 ug
Disease State	Cancer
Diagnosis Category	Colorectal Cancer
Tissue	Colon
Frozen Sample ID	FR00018DBE
Case ID	CU0000005634
Prep	Not applic
Vial ID	FR00018DBE
Tissue of Origin	Colon: right
Site of Finding	Colon: right
Appearance	Tumor
Sample Diagnosis (from Pathology Verification)	Adenocarcinoma of colon
Normal	10%
Lesional	0%
Tumor	60%
Tumor Hypo/Acellular Stroma	0%
Tumor Hypercellular Stroma	20%
Necrosis	10%
Pathology Verification notes from H&E review	Tumor Stroma (Cellular): Inflammatory cells
CASE Diagnosis (from medical center path report)	Adenocarcinoma of colon
Age	61
Gender	Female
Tumor Grade	AJCC G1: Well differentiated
Minimum Stage Grouping	III
T (Extent of Primary Tumor)	T4
N (Lymph node metastasis)	N1
M (Distant metastasis)	MX
Storage & Handling	The lysate is shipped with dry ice. Store the tissue lysate at -80°C. These products are stable for one year from date of shipping when stored at -80°C without repeated freeze/thaws.

Submit a product review to Biocompare.com





Submit a review of this product to Biocompare.com to receive a $20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.

BOSTER BIOLOGICAL TECHNOLOGY
3942 B Valley Ave, Pleasanton, CA 94566
888-466-3604 | support@bosterbio.com | www.bosterbio.com

For Research Use Only. Not for use in diagnostic procedures.