BOSTER BIOLOGICAL TECHNOLOGY
3942 B Valley Ave, Pleasanton, CA 94566
888-466-3604 | support@bosterbio.com | www.bosterbio.com

Anti-ERBB1 / EGFR / HER1 Reference Antibody (panitumumab)

Catalog Number: M00023-11

Submit a product review to Biocompare.com





Submit a review of this product to Biocompare.com to receive a $20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Overview


Product Name	Anti-ERBB1 / EGFR / HER1 Reference Antibody (panitumumab)
Reactive Species	Human
Description	Boster Bio Anti-ERBB1 / EGFR / HER1 Reference Antibody (panitumumab) (Catalog # M00023-11). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Human. Endotoxin: < 0.059EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Formulation	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Storage Instructions	Maintain refrigerated at 2–8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	P00533

Technical Details


Isotype	IgG2SA
Form	Liquid

BOSTER BIOLOGICAL TECHNOLOGY
3942 B Valley Ave, Pleasanton, CA 94566
888-466-3604 | support@bosterbio.com | www.bosterbio.com


Anti-ERBB1 / EGFR / HER1 Reference Antibody (panitumumab) (M00023-11) Images




Anti-ERBB1/EGFR/HER1 Reference Antibody (panitumumab) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%



The purity of Anti-ERBB1/EGFR/HER1 Reference Antibody (panitumumab)is more than 98.8%



Immobilized human EGFR His at 2 &mug/mL can bind Anti-ERBB1/EGFR/HER1 Reference Antibody (panitumumab)



NCI-H292 cells were stained with Anti-ERBB1/EGFR/HER1 Reference Antibody (panitumumab) and negative control protein respectively

For Research Use Only. Not for use in diagnostic procedures.