BOSTER BIOLOGICAL TECHNOLOGY
3942 B Valley Ave, Pleasanton, CA 94566
888-466-3604 | support@bosterbio.com | www.bosterbio.com

TCTN1 (NM_001082537) Human Recombinant Protein

Catalog Number: PROTQ2MV58

Submit a product review to Biocompare.com





Submit a review of this product to Biocompare.com to receive a $20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Overview

Product Name	TCTN1 (NM_001082537) Human Recombinant Protein
Description	Recombinant protein of human tectonic family member 1 (TCTN1), transcript variant 2
Size	20 μg
Tag	C-Myc/DDK
Form	Frozen Solution in PBS Buffer
Source	HEK293T
Formulation	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol

Concentration

>50 ug/mL as determined by microplate BCA method

Storage

Store at -80°C. Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles. (Ship on dry ice.)

Purity

> 80% as determined by SDS-PAGE and Coomassie blue staining

Amino Acid Sequence

MRPRGLPPLLVVLLGCWASVSAQTDATPAVTTEGLNSTEAALATFGTFPSTRPPGTPRAPGPSSGPRPTPVTDVAVLCVCDLSPAQCDINCCCDPDCSSV
DFSVFSACSVPVVTGDSQFCSQKAVIYSLNFTANPPQRVFELVDQINPSIFCIHITNYKPALSFINPEVPDENNFDTLMKTSDGFTLNAESYVSFTTKLDIPT
AAKYEYGVPLQTSDSFLRFPSSLTSSLCTDNNPAAFLVNQAVKCTRKINLEQCEEIEALSMAFYSSPEILRVPDSRKKVPITVQSIVIQSLNKTLTRREDTDVL
QPTLVNAGHFSLCVNVVLEVKYSLTYTDAGEVTKADLSFVLGTVSSVVVPLQQKFEIHFLQENTQPVPLSGNPGYVVGLPLAAGFQPHKGSGIIQTTNRY
GQLTILHSTTEQDCLALEGVRTPVLFGYTMQSGCKLRLTGALPCQLVAQKVKSLLWGQGFPDYVAPFGNSQAQDMLDWVPIHFITQSFNRKDSCQLPG
ALVIEVKWTKYGSLLNPQAKIVNVTANLISSSFPEANSGNERTILISTAVTFVDVSAPAEAGFRAPPAINARLPFNFFFPFV

BOSTER BIOLOGICAL TECHNOLOGY
3942 B Valley Ave, Pleasanton, CA 94566
888-466-3604 | support@bosterbio.com | www.bosterbio.com

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

TCTN1 (NM_001082537) Human Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.