BOSTER BIOLOGICAL TECHNOLOGY
3942 B Valley Ave, Pleasanton, CA 94566
888-466-3604 | support@bosterbio.com | www.bosterbio.com

TUG (ASPSCR1) (NM_024083) Human Recombinant Protein

Catalog Number: PROTQ9BZE9

Submit a product review to Biocompare.com





Submit a review of this product to Biocompare.com to receive a $20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Overview

Product Name	TUG (ASPSCR1) (NM_024083) Human Recombinant Protein
Description	Recombinant protein of human alveolar soft part sarcoma chromosome region, candidate 1 (ASPSCR1)
Size	20 μg
Tag	C-Myc/DDK
Form	Frozen Solution in PBS Buffer
Source	HEK293T
Formulation	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol

Concentration

>50 ug/mL as determined by microplate BCA method

Storage

Store at -80°C. Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles. (Ship on dry ice.)

Purity

> 80% as determined by SDS-PAGE and Coomassie blue staining

Amino Acid Sequence

MAAPAGGGGSAVSVLAPNGRRHTVKVTPSTVLLQVLEDTCRRQDFNPCEYDLKFQRSVLDLSLQWRFANLPNNAKLEMVPASRSREGPENMVRIALQL
DDGSRLQDSFCSGQTLWELLSHFPQIRECLQHPGGATPVCVYTRDEVTGEAALRGTTLQSLGLTGGSATIRFVMKCYDPVGKTPGSLGSSASAGQAAAS
APLPLESGELSRGDLSRPEDADTSGPCCEHTQEKQSTRAPAAAPFVPFSGGGQRQGGPPGPTRPLTSSSAKLPKSLSSPGGPSKPKKSKSGQDPQQEQE
QERERDPQQEQERERPVDREPVDREPVVCHPDLEERLQAWPAELPDEFFELTVDDVRRRLAQLKSERKRLEEAPLVTKAFREAQIKEKLERYPKVALRVL
FPDRYVLQGFFRPSETVGDLRDFVRSHLGNPELSFYLFITPPKTVLDDHTQTLFQANLFPAALVHLGAEEPAGVYLEPGLLEHAISPSAADVLVARYMSRA
AGSPSPLPAPDPAPKSEPAAEEGALVPPEPIPGTAQPVKRSLGKVPKWLKLPASKR

BOSTER BIOLOGICAL TECHNOLOGY
3942 B Valley Ave, Pleasanton, CA 94566
888-466-3604 | support@bosterbio.com | www.bosterbio.com

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

TUG (ASPSCR1) (NM_024083) Human Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.